



ORTHOTICS

Collection



Welcome to Our Orthotic Catalog

03-04

Who develops our AFOs?

Trulife's Carbon Fiber Dynamic AFOs

05-11

Ankle Foot Orthotics

Pace Carbon Fiber AFOs

Matrix Ankle Foot Orthotics

12-13

Cervical

Lerman Minerva

Lerman Non-Invasive Halo

14-18

Spinal

Universal CASH Hyperextension

T34 Hyperextension Orthosis

T39 Hyperextension Orthosis

FAST-WRAP

19-19

Maternity Supports

7205 Embrace Moderate Support Maternity Belt

7210 Embrace Ultimate Support Maternity Belt

22-22

Knee Orthoses

Lenox Hill Hinged Swedish Knee Cage

Lenox Hill Swedish Knee Cage

21-22

Pediatric

Pectus Carinatum Orthosis

Pectus Pro Pectus Carinatum Orthosis

Pectus Pro-S Pectus Carinatum Orthosis

23-23

Contact Us



Who develops our AFOs?

Meet Abraham

Trulife Research and Development
Engineer

Abraham, originally from Mexico, began working for Trulife in 2020. With more than a decade of professional experience, he has successfully led projects for the aerospace, automotive, and medical device industries.

Abraham leads from concept through validation, manufacturing, and global launch of the development of Class I medical devices. He is responsible for the design and development of some of our acclaimed Orthotic devices, including the Pace Line, the Matrix range and the Pectus Pro Chest Braces.



Trulife's Carbon Fiber Dynamic AFOs

The Matrix Family has been developed to provide a dynamic response for patients with drop foot, limb proprioception deficit, mild proximal deficit as well as for partial foot amputations.

The Pace Line combines advanced carbon fiber engineering with refined fit and function.

Designed to enhance gait efficiency, reduce fatigue, and deliver long-lasting support for patients with foot and ankle deficits.



Model	FH400	FH450	FH500	FH200	FH250	FH260	FH270	FH225
Name	M-Pace	L-Pace	A-Pace	Matrix	Matrix Max	Matrix Max 2 (Trimmable)	Matrix SuperMax	Matrix Max Raw
Sizes	S, M, L, XL	S, M, L, XL	S/M, L/XL	S, M, L	S, M, L, XL	S, M, L, XL	S, M, L, XL	ONE SIZE
Anterior/Posterior	Posterior	Posterior	Anterior	Anterior	Anterior	Anterior	Anterior	Anterior
Product Image								

Trulife AFO Rigidity Chart



Pace Carbon Fiber AFOs

A-Pace Anterior Lateral Strut AFO

Anterior-lateral strut AFO offering intermediate stiffness for an optimal balance of support, flexibility and energy return - ideal for an active lifestyle.

- Fewer joints to enhance strength and reliability.
- Magnetic closures allow for easy, one-handed operation. Ideal for users with reduced hand mobility or dexterity.
- Ergonomically designed footplate with curved heel for a more anatomical fit.
- A narrower footplate allows compatibility with a wider variety of footwear.
- Cushioned soft padding that is stylish and breathable. Providing all-day comfort.

Indications:

Low to high activity patients.
Mild to moderate footdrop.
With or without mild to moderate foot spasticity.
With weak muscles or with loss of calf muscle.
For unilateral and bilateral patients.



Customizable height with a flexible, high-strength front guard that conforms comfortably to the calf.



Made from layered Formula 1 and aerospace-grade carbon fiber.



Fewer SKUs with two size options with a trim-to-fit footplate.

Selection

Model No. Right	Model No. Left	Size	Height	Foot Length
FH500R01	FH500L01	Small/Medium	36.3-39.7cm / 14.3-15.6in	23-27cm / 9.1-10.6in
FH500R02	FH500L02	Large/X-Large	37.3-40.7cm / 14.7-16in	26-30cm / 10.2-11.8in

Specifications

A-Pace	
Material	Carbon Fiber
Warranty	Hardware: 2 years Soft goods: 90 days

SCAN ME! >>>



Pace Carbon Fiber AFOs

L-Pace Posterior Lateral Strut AFO

- Posterior lateral strut of intermediate stiffness—providing an optimal mix of support, flexibility, and energy return.
- Constructed from aerospace-grade carbon fiber materials for lightweight strength and durability.
- Ergonomic footplate with a curved heel for enhanced fit and comfort.
- A narrower footplate increases compatibility with a variety of footwear.
- Reverse-pull, hook-and-loop strap for a secure and comfortable fit.
- Flexible, high-strength posterior band that conforms easily to different calf shapes.

Indications:

Low to high activity patients with unilateral or bilateral drop foot



Ergonomic, slim footplate is designed to accommodate a wider variety of footwear.



Breathable padded cuff with reverse-pull strap.



Aerospace-grade carbon fiber materials for strength without bulk.

Selection

Model No. Right	Model No. Left	Size	Height	Foot Length	Shoe Size US	Shoe Size EU	Shoe Size UK
FH450R02	FH450L02	Small	36.5cm / 14.4in	<23.6cm / <9.3in	5/6	36/38	4/5.5
FH450R03	FH450L03	Medium	36.5cm / 14.4in	25cm / 9.8in	7/8	39/40	6/7
FH450R04	FH450L04	Large	38.5cm / 15.2in	26.4cm / 10.4in	9/10	41/43	7.5/9.5
FH450R05	FH450L05	X-Large	38.5cm / 15.2in	27.4cm / 10.8in	10>	43>	9.5>

Specifications

L-Pace	
Material	Carbon Fiber
Warranty	Hardware: 2 years Soft goods: 90 days

SCAN ME! >>>



Pace Carbon Fiber AFOs

M-Pace Posterior Medial Strut AFO

- Posterior medial strut of intermediate stiffness —providing an optimal mix of support, flexibility, and energy return.
- Constructed from aerospace-grade carbon fiber materials for lightweight strength and durability.
- Ergonomic footplate with a curved heel for enhanced fit and comfort.
- A narrower footplate increases compatibility with a variety of footwear.
- Reverse-pull, hook-and-loop strap for a secure and comfortable fit.
- Flexible, high-strength posterior band that conforms easily to different calf shapes.

Indications:

Low to high activity patients with unilateral or bilateral drop foot



Ergonomic, slim footplate is designed to accommodate a wider variety of footwear.



Breathable padded cuff with reverse-pull strap.



Aerospace-grade carbon fiber materials for strength without bulk.

Selection

Model No. Right	Model No. Left	Size	Height	Foot Length	Shoe Size US	Shoe Size EU	Shoe Size UK
FH400R02	FH400L02	Small	36.5cm / 14.4in	<23.6cm / <9.3in	5/6	36/38	4/5.5
FH400R03	FH400L03	Medium	36.5cm / 14.4in	25cm / 9.8in	7/8	39/40	6/7
FH400R04	FH400L04	Large	38.5cm / 15.2in	26.4cm / 10.4in	9/10	41/43	7.5/9.5
FH400R05	FH400L05	X-Large	38.5cm / 15.2in	27.4cm / 10.8in	10>	43>	9.5>

Specifications

M-Pace	
Material	Carbon Fiber
Warranty	Hardware: 2 years Soft goods: 90 days

SCAN ME! >>>



Pace Carbon Fiber AFOs

The Pace Line Replacement Pad & Strap Kits

A-Pace Kit

Replacement Pad & Strap Kit for A-Pace - Anterior Lateral Strut AFO



Specifications		
A-Pace Replacement Pad & Strap Kit		
Item	Description	Quantity
1	Pad	1
2	Straps & Magnetic Buckles	2
3	Self-Adhesive Hooks	10

L-Pace & M-Pace Kit

Replacement Pad & Strap Kit suitable for both the L-Pace & M-Pace.



Specifications		
L-Pace & M-Pace Replacement Pad & Strap Kit		
Item	Description	Quantity
1	450mm length velour BTB strap	1
2	50mm tab with Trulife logo	1
3	Back to Back 48mm (Hook)	1
4	Pace Pad	1
5	70mm length black micro 3 hook adhesive strips	4
6	25mm hook and loop dots	4

Matrix Ankle Foot Orthotics

Matrix

- Energy storing, composite footplate and strut.
- Lateral strut is anterior to malleolus.
- Open calcaneal design.
- Custom-fit, anterior shell and footplate.

Indications:

Low activity patients with simple, unilateral drop foot, maximum weight of 113 kg / 250 lb



Selection				
Model No. Right	Model No. Left	Size	Height	Foot Length
FH200R02	FH200L02	Small	39-42cm / 15.5-16.5in	<25cm / <9.8in
FH200R03	FH200L03	Medium	41-43cm / 16-17in	24-27cm / 9.4-10.6in
FH200R04	FH200L04	Large	43-46cm / 17-18in	25-29cm / 9.8-11.4in

Specifications	
Matrix	
Material	Carbon Composite
Includes	Hex Key, Tibial Pad x2
Warranty	Hardware: 2 years Soft goods: 90 days

Matrix Max

- Lightweight, carbon composite frame with enhanced stability.
- Reinforced lateral strut that is anterior to malleolus.
- Open calcaneal design.
- Custom-fit, anterior shell and footplate.

Indications:

Low to high activity patients with unilateral or bilateral drop foot; suitable for partial foot amputations.



Selection				
Model No. Right	Model No. Left	Size	Height	Foot Length
FH250R02	FH250L02	Small	39-42cm / 15.5-16.5in	<25cm / <9.8in
FH250R03	FH250L03	Medium	41-43cm / 16-17in	24-27cm / 9.4-10.6in
FH250R04	FH250L04	Large	43-46cm / 17-18in	25-29cm / 9.8-11.4in
FH250R05	FH250L05	X-Large	43-46cm / 17-18in	26-31cm / 10.2-12.2in

Specifications	
Matrix Max	
Material	Carbon Composite
Includes	Hex Key, Tibial Pad x2
Warranty	Hardware: 2 years Soft goods: 90 days

Matrix Ankle Foot Orthotics

Matrix Max 2 (Trimmable)

- Easy-Trim footplate allows fitting of up to 3 shoe sizes with one product.
- Advanced reinforced lateral strut for greater stability, durability and energy return.
- Narrower foot plate in the medial arch, allowing for an easier shoe fit.
- Heat-mouldable, trimmable and height-adjustable anterior shell.

Indications:

Low to high activity patients with unilateral or bilateral drop foot; suitable for partial foot amputations.



Selection				
Model No. Right	Model No. Left	Size	Height	Foot Length
FH260R02	FH260L02	Small	39-42cm / 15.5-16.5in	22-24cm / 8.5-9.5in
FH260R03	FH260L03	Medium	41-43cm / 16-17in	24-26cm / 9.5-10.2in
FH260R04	FH260L04	Large	43-46cm / 17-18in	26-29cm / 10.2-11.5in
FH260R05	FH260L05	X-Large	43-46cm / 17-18in	29-31cm / 11.5-12.2in

Specifications	
Matrix Max 2 (Trimmable)	
Material	Carbon Composite
Width at Medial Arch	Small: 6cm / 2.3in Medium: 6.3cm / 2.5in Large: 6.6cm / 2.6in X-Large: 6.9cm / 2.7in
Includes	Hex Key, Tibial Pad x2
Warranty	Hardware: 2 years Soft goods: 90 days

Matrix SuperMax

- Enhanced integration of the foot plate and strut for added durability and support.
- Footplate is uniquely designed to allow greater energy reflection at toe off.
- Streamlined to fit in shoes better.
- Shank angle increased to 82° to improve knee extension movement at mid-stance.
- Heat-mouldable, trimmable and height-adjustable anterior shell.

Indications:

Low to high activity patients with unilateral or bilateral drop foot, peroneal nerve palsy, incomplete paralysis, CVA and CMT.



Selection				
Model No. Right	Model No. Left	Size	Height	Foot Length
FH270R02	FH270L02	Small	39-42cm / 15.5-16.5in	<25cm / <9.8in
FH270R03	FH270L03	Medium	41-43cm / 16-17in	24-27cm / 9.4-10.6in
FH270R04	FH270L04	Large	43-46cm / 17-18in	25-29cm / 9.8-11.4in
FH270R05	FH270L05	X-Large	43-46cm / 17-18in	26-31cm / 10.2-12.2in

Specifications	
Matrix SuperMax	
Material	Carbon Composite
Includes	Hex Key, Tibial Pad x2
Warranty	Hardware: 2 years Soft goods: 90 days

Matrix Ankle Foot Orthotics

Matrix Max Raw

- Lightweight, carbon composite.
- Ready for practitioner to create a custom anterior AFO, or partial foot prosthesis.
- One size.

Indications:

Drop foot, CMT, MS, ALS, off-loading PTB and any other creative ways that provide optimal outcomes for patients.



Selection				
Model No. Right	Model No. Left	Size	Height	Foot Length
FH225R	FH225L	One Size	37cm / 14.5in	29cm / 11.5in

Specifications	
Matrix Max Raw	
Material	Carbon Composite
Includes	White Shrink Tube, Stainless Steel Inserts, M4 Screws, Black Shrink Tube x2
Warranty	Hardware: 2 years

Matrix Composite Struts

- Prepreg composite struts for manufacturing custom posterior leaf spring type AFOs.
- Available in 2 widths, each with two levels of flexibility (Medium, Firm)
- Length of 18" (45 cm), which is to be trimmed to the required length.
- Provides excellent mediolateral stability similar to a rigid ankle type AFO while still allowing graded plantarflexion and dorsiflexion.

Indications:

Patients requiring a custom posterior leaf spring AFO for a variety of conditions.



Selection		
Model No.	Flexibility	Width
FH255-M	Medium	30mm
FH265-M	Medium	40mm
FH255-F	Firm	30mm
FH265-F	Firm	40mm

Specifications	
Matrix Composite Struts	
Material	Carbon Composite
Length	45cm / 18in
Warranty	1 Year

Lerman Minerva

- Manufactured with non-ferrous material.
- Used to support and Stabilize cervical spine and upper thoracic spine after injury or surgery.
- Replacement pad sets available, (Includes: Back Pad, Chin Pad (2), Chest Pad, Occipital Pad, Head Band Assembly)

Indications:

C2-T3 stable fractures, muscular/ligamentous injury, degenerative joint disease, cervical pain/nerve impingement.

Contraindications:

Unstable fractures, type II dens fractures.



Selection

Model No.	Product Description	Size	Height
A193PED	Lerman Minerva	Paediatric	23-28 / 9-11in
A193SM	Lerman Minerva	Small	28-33cm / 11-13in
A193REG	Lerman Minerva	Regular	33-41cm / 13-16in
O11394	Replacement Pad Set	Paediatric	n/a
O10902	Replacement Pad Set	Small	n/a
O10901	Replacement Pad Set	Regular	n/a

Specifications

Lerman Minerva	
Material	Polyethylene / Aluminium / Latex Free
Includes	2 Pad Sets
Warranty	90 Days

Cervical

Lerman Non-Invasive Halo

- Secure cervical stabilization.
- Colour-coded buckles for ease of application.
- Designed to be donned/doffed in supine position to simplify fit.
- Removable anterior plate for emergency Cardiopulmonary Resuscitation.
- Polyethylene / composite plate material.
- Manufactured with non-ferrous materials.
- One size fits most.

Indications:

Cervical spine immobilization of stable, non-displaced or minimally displaced fractures.



Selection	
Model No.	Size
A37802	One Size Fits Most

Specifications	
Lerman Non-Invasive Halo	
Material	Polyethylene / Composite Plate / Latex Free
Includes	Pad Sets
Warranty	90 Days

Universal CASH Hyperextension

- Universal posterior strap can be cut to length from 76 cm/30in up to 122cm/ 48in.
- Universal posterior pad can be oriented vertically or horizontally.
- Fully adjustable height and width.
- Easy application and removal.
- Closed foam padding for easy cleaning .

Indications:

Stable compression fractures of T7-L2, osteoporosis, kyphosis and osteoarthritis.



Selection		
Model No.	Product Description	Size
CASH100	Fixed Sternal and Pubic Pad	One Size Fits Most
CASH300S	Hinged Sternal and Pubic Pad	One Size Fits Most
CASH300P	Hinged Pectoral and Pubic Pad	One Size Fits Most

Specifications	
Universal CASH Hyperextension	
Material	Aluminium / Low Density Polyethylene / Latex Free
Warranty	90 Days

Universal CASH Hyperextension Accessories/Replacement Parts

Selection	
Model No.	Product Description
CASH103	CASH Long Horizontal Bar
CASH105	CASH Long Vertical Bar
CASH106	CASH Short Horizontal Bar
CASH107	CASH Short Vertical Bar
CASH116	CASH Posterior Strap 48in
CASH125	CASH Kodel Sternal Pad Cover
CASH126	CASH Kodel Pectoral Pad Cover (Pair)
CASH135	CASH Front Closure Kit (No Strap)



Spinal

T34 Hyperextension Orthosis with Pubic Pad

- Anatomically designed pubic pad and pelvic bar avoids pressure and prevents migration.
- Sternal pad swivels to eliminate pressure.
- Latch-lock clamp closure allows easy donning and doffing.
- Adjustable height and width.

Indications:

Vertebral compression fractures, osteoporosis, vertebral arthritis and post-surgical support.



Selection		
Model No.	Size	Hip Circumference
T34-02	Small	76-83cm / 30-32.5in
T34-03	Medium	83-94cm / 32.5-37in
T34-04	Large	94-104cm / 37-41in
T34-05	X-Large	104-114cm / 41-45in

Specifications	
T34 Hyperextension Orthosis with Pubic Pad	
Material	Aluminium / Latex Free
Brace Height	Small: 43-47cm / 17-18.5in Medium: 45-50cm / 17.5-20in Large: 48-54cm / 19-21.5in X-Large: 53-59cm / 21-23in
Warranty	90 Days

T39 Hyperextension Orthosis with Articulating, Adjustable Pelvic Band

- Articulating, adjustable pelvic band.
- Posterior pad may be used horizontally or vertically.
- Swivel sternal pad for patient comfort.
- Latch-lock clamp closure allows ease of donning and doffing.
- Available in regular (T39) or short (T39S)
- M/L and height adjustable.
- Removable, washable pads .

Indications:

Osteoporosis, thoracolumbar fractures T6-L1 and post surgical support.



Selection				
Model No.	Product Description	Size	Hip Circumference	Pelvic Band Length
T39-02	T39 Regular	Small / Medium	76-94cm / 30-37in	43-48cm / 17-19in
T39-04	T39 Regular	Large / X-Large	94-114cm / 37-45in	53-59cm / 21-23.5in
T39S-02	T39 Short	Small / Medium	76-94cm / 30-37in	43-48cm / 17-19in
T39S-04	T39 Short	Large / X-Large	94-114cm / 37-45in	53-59cm / 21-23.5in

Specifications	
T39 Hyperextension Orthosis with Articulating, Adjustable Pelvic Band	
Material	Aluminium / Latex Free
Brace Height	(Regular) Small/Medium: 43-50cm / 17-19.5in Large/X-Large: 50-60cm / 19.5-23.5in (Short) Small/Medium: 39-43cm / 15-17in Large/X-Large: 43-50cm / 17-19.5in
Warranty	90 Days

Spinal

FAST-WRAP Thoracolumbar Support

- Compression and stabilization of thoracic, lumbar, sacral and abdominal regions.
- Lower cinch straps adjust level of support.
- Padded, adjustable shoulder straps.
- Contoured, posterior 43 cm / 17 in steel stays provide additional support.

Indications:

Back sprains/strains, osteoarthritis, mild osteoporosis, posture control.



Selection		
Model No.	Size	Hip Circumference
0078002	Small	81-91cm / 32-36in
0078003	Medium	97-107cm / 38-42in
0078004	Large	112-122cm / 44-48in
0078005	X-Large	127-132cm / 50-52in
0078006	XX-Large	137-142cm / 54-56in

Specifications	
FAST-WRAP Thoracolumbar Support	
Material	Nylon / Polyurethane
Stays	Contoured Steel
Back Height	53cm / 21in
Front Height	18cm / 7in
Warranty	90 Days

FAST-WRAP Semi-Pendulous Support with Spring Steel Stays

- Compresses and Stabilizes lumbar, sacral and abdominal regions.
- Cupped front accommodates semi-pendulous abdomen.
- Steel stays offer additional support.

Indications:

Low back pain; ideal for geriatric or post-natal patients.



Selection		
Model No.	Size	Hip Circumference
0035702	Small	86-91cm / 34-36in
0035703	Medium	97-102cm / 38-40in
0035704	Large	107-112cm / 42-44in
0035705	X-Large	117-122cm / 46-48in
0035706	XX-Large	127-132cm / 50-52in

Specifications	
FAST-WRAP Semi-Pendulous Support with Spring Steel Stays	
Material	Nylon / Polyurethane / Latex Free
Stays	Steel
Back Height	25cm / 10in
Front Height	23cm / 9in
Warranty	90 Days

PT132 Platinum Lumbosacral Support

- Adjustable, discreet lumbosacral stabilization and abdominal support.
- Antibacterial Trutex fabric promotes freshness.
- Breathable stretch foam.

Indications:

Low back sprains/strains, mild lumbar disc injury.



Selection		
Model No.	Size	Hip Circumference
PT13202	Small	74-84cm / 29-33in
PT13203	Medium	85-100cm / 34-39in
PT13204	Large	101-115cm / 40-45in
PT13205	X-Large	116-129cm / 46-51in

Specifications	
PT132 Platinum Lumbosacral Support	
Material	Polyester / Polyurethane / Spandex / Latex Free
Posterior Stays	Plastic x4 (24cm / 9.5in)
Back Height	28cm / 11in
Front Height	19cm / 7.6in
Warranty	90 Days

PT137 Platinum Thoracolumbar Orthosis

- Adjustable, discreet thoracic, lumbar, sacral and abdominal support.
- Stays and boning provide additional back support.
- Antibacterial Trutex fabric promotes freshness.
- Breathable stretch foam.

Indications:

Back sprains/strains, osteoarthritis, mild osteoporosis, posture control.



Selection		
Model No.	Size	Hip Circumference
PT13702	Small	75-85cm / 29.5-33.4in
PT13703	Medium	85-95cm / 33.4-37.4in
PT13704	Large	95-105.5cm / 37.5-41.5in
PT13705	X-Large	105.5-115.5cm / 41.5-45.5in
PT13706	XX-Large	115.5-126cm / 45.5-49.5in

Specifications	
PT137 Platinum Thoracolumbar Orthosis	
Material	Polyester / Polyurethane / Spandex
Stays	Plastic x4, 38cm / 15in (x2), 48cm / 19in (x2)
Back Height	56cm / 22in
Front Height	23cm / 9in
Warranty	90 Days

Maternity Supports

7205 Embrace Moderate Support Maternity Belt

- Moderate abdominal and lumbar support.
- Evenly transfers the weight of the abdomen to the spine.

Indications:

Back pain associated with pregnancy.



Selection		
Model No.	Size	Hip Circumference
0720502	Small	86-102cm / 34-40in
0720503	Medium	97-112cm / 38-44in
0720504	Large	107-122cm / 42-48in
0720505	X-Large	117-132cm / 46-52in
0720506	XX-Large	127-142cm / 50-56in

Specifications	
7205 Embrace Moderate Support Maternity Belt	
Material	Cotton / Rubber / Polyester / Latex Free
Back Height	15cm / 6in
Front Height	9cm / 3.5in
Warranty	90 Days

7210 Embrace Ultimate Support Maternity Belt

- Maximum abdominal and lumbar support.
- Evenly transfers the weight of the abdomen to the spine.
- Contours to expanding abdomen.
- Low temperature, heat-mouldable, rigid back support panel available.

Indications:

Back pain associated with pregnancy.



Heat-mouldable
Rigid Back
Support Panel
Available ▼

Selection		
Model No.	Size	Hip Circumference
0721002	Small	86-102cm / 34-40in
0721003	Medium	97-112cm / 38-44in
0721004	Large	107-122cm / 42-48in
0721005	X-Large	117-132cm / 46-52in

Specifications	
7210 Embrace Ultimate Support Maternity Belt	
Material	Cotton / Rubber / Polyester / Latex Free
Back Height	20cm / 8in
Front Height	15cm / 6in
Warranty	90 Days

Lenox Hill Hinged Swedish Knee Cage

- Lightweight, aluminium frame with padded straps.
- Prevents hyperextension of knee joint.
- Will not protrude above or below knee while seated .

Indications:

Knee hyperextension, genu recurvatum.



Selection		
Model No.	Size	Knee Circumference*
A24700	Small	28-38cm / 11-15in
A24701	Regular	38-53cm / 15-21in

*Measure 10cm / 4in above knee center

Specifications	
Lenox Hill Hinged Swedish Knee Cage	
Material	Aluminium (Frame/Joint)
Joint Type	Hinged
Thigh/Calf Band	Anterior/Posterior
Warranty	90 Days

Lenox Hill Swedish Knee Cage

- Lightweight, vinyl coated, aluminium frame.
- Moulds easily.
- Adjustable, posterior pressure pad.
- Permits knee flexion, restricts extension.

Indications:

Knee hyperextension, genu recurvatum.



Selection		
Model No.	Size	Knee Circumference*
0866802	Small	25-30cm / 10-12in
0866803	Medium	33-35cm / 13-14in
0866804	Large	38-41cm / 15-16in
0866805	X-Large	42-46cm / 17-18in

*Measure 10cm / 4in above knee center

Specifications	
Lenox Hill Swedish Knee Cage	
Material	Powder Coated Aluminium
Thigh/Calf Band	Anterior/Anterior
Warranty	90 Days

Pectus Carinatum Orthosis

- Provides increased anterior/posterior directed corrective compression forces.
- Bendable metal components.
- Washable foam pads increase patient comfort.
- Ratchet closure.
- One size fits most.

Indications:

Abnormal chest shape, protrusion of both the sternum (breast bone) and the cartilage portion of the ribs, bowed chest (pigeon breast)



Selection		
Model No.	Size	Chest Circumference
PECTUS	One Size Fits Most	66-89cm / 26-35in

Specifications	
Pectus	
Material	Low Density Polyethylene/ Aluminium/ Padding / Latex Free
Warranty	90 Days



SCAN TO VIEW FITTING VIDEO

Pectus Pro Pectus Carinatum Orthosis

- Provides increased anterior/posterior directed corrective compression forces.
- New and improved design.
- Wider size range (29-41") to accommodate larger chests.
- New ratchet closure for greater durability.
- Optional, soft, shoulder straps included for additional stabilization during wear.
- 2.5" of lateral adjustment for asymmetric Pectus Carinatum.
- Removable, washable and breathable pads for increased comfort.
- Slide adjustment bars and screws.
- One size fits most.

Indications:

Abnormal chest shape, protrusion of both the sternum (breastbone) and the cartilage portion of the ribs, and bowed chest (Pigeon chest)



Selection		
Model No.	Size	Chest Circumference
PECTUS PRO	One Size Fits Most	74-104cm / 29-41in

Specifications	
Pectus Pro	
Material	Low Density Polyethylene/ Aluminium/ Padding / Latex Free
Warranty	90 Days

Pectus Pro-S Pectus Carinatum Orthosis

The new Pectus Pro-S offers the same effective compression forces for Pectus Carinatum and other chest anomalies as our original orthosis, now redesigned for smaller patients.

- Shorter adjustment bars, a smaller back panel, and improved comfort.
- Smaller chest size range of 63-102cm.
- Direct anterior-posterior corrective compression forces.
- Durable ratchet closure.
- 2.5" of lateral adjustment for asymmetric Pectus Carinatum.
- Black anodised slide adjustment bars and screws.
- Optional, soft, shoulder straps included.
- Removable, breathable, washable pads.

Indications:

Abnormal chest shape, protrusion of both the sternum (breastbone) and the cartilage portion of the ribs, and bowed chest (Pigeon chest)



Selection		
Model No.	Size	Chest Circumference
PECTUS PRO-S	~ 63 cm	~ 102 cm

Specifications	
Pectus Pro-S	
Warranty	90 Days

Any questions?

Contact US

Whether you have a question about our products, need assistance, or simply want to share your feedback, our team is ready to assist you.

Scan the QR code below for contact details.





Pace
CARBON FIBER AFOS