

Click Medical products are patent protected.

For full list see www.clickmedical.co/patents

RevoSurface®

SKU	Product	Description	Picture
RS-100-07-01	RevoSurface® Tool Kit with Click® Reel	This single RevoSurface Tool Kit with the Click Reel allows you to test adjustable designs and fit of check sockets, 3D printed devices, and on orthoses. It can also be used to solve fit issues on existing rigid devices. The Click Reel innovation creates adjustable devices with a lower streamlined profile, micro/macro adjustment that incrementally tightens and loosens with increased durability. Kit includes: 1x Click Reel, 2x Surface Collar, 2m HC Lace, 1.8m Lamination Tubing, 1x Metal Lace Feeder, 3x Plastic Lace Feeders, 1x Reel Tool, 6x Surface Guides, 3x Releasing Guides, 4x Chicago Screw assemblies User Weight Limit: 300lb (135kg)	
RP-001	HC Lace Kit	2m of HC (High Capacity) Lace and Click Reel components for RevoSurface Tool Kit, RevoLock Align, and RevoLock 4-Hole Kits with the Click Reel.	

RP-000-07	Click® Reel and HC Lace Kit	Replacement Click Reel and 2m of HC (High Capacity) Lace for RevoLock Align, RevoLock 4-Hole and RevoSurface Tool Kits only. This is not for replacing a BOA® dial. Please see RP0100-320-05 for replacing a BOA® dial.	 The image shows the components of the Click Reel and HC Lace Kit: a black Click Reel, a coiled black HC (High Capacity) Lace, a small grey plastic piece labeled 'CLICK', and a coiled green HC Lace.
RP-101	Surface Guide 6-Pack	The 6-Pack of Surface Guides allows a practice to have on hand additional guides for easy replacement and configuration changes. For RevoSurface Tool Kits with Click Reel only.	 The image shows six black Surface Guides, which are small, circular, and have a central hole.
RP-102	Surface Collar	2x Click Reel mounting collar and 4x Chicago Screw assemblies. For surface mount only. Not for use with lamination.	 The image shows two black Click Reel mounting collars, which are circular with four mounting holes, and four black Chicago Screw assemblies.
RP-105	Surface Releasing Tab 3-Pack	3-Pack Surface Releasing Tabs for converting a device from a fixed lace path to a releasing lace path to fully open the device or to replace existing Surface Releasing Tabs.	 The image shows three black Surface Releasing Tabs, which are small, L-shaped plastic pieces.
RP-201	Lamination Tube - 1.8m	Replacement Lamination Tube (1.8m) for all Revo Kits with Click Reel and BOA® dials.	 The image shows a coiled white Lamination Tube, which is a flexible, translucent plastic tube.

BL-001	HC Lace - 25m Spool	25m spool of HC (High Capacity) Lace assures you will always have replacement lace on hand. Cost effective, sustainable and less waste. Used in the RevoSurface, RevoLock Align, RevoLock 4-Hole, and RevoLock Lanyard V2 & V3 with BOA Dial.	
AC-100-05	Comfort Grip Ring - Black 5-Pack	Provide a soft-touch feel and 10% increased turning power to the existing Click Reel. Great for patients with dexterity issues.	
AC-101-05	Power Grip Ring - Black 5-Pack	Provide a soft-touch feel and 20% increased turning power to the existing Click Reel. Great for patients with dexterity and hand strength issues.	
AC-000-00	Click Care Kit™	The Click Care Kit enables you to service your clients in the office or on the road with everything readily accessible for maintaining, repairing or replacing Click components. Contents: 1x Surface Collar Kit, 2x HD Lace Kits, 2x HC Lace Kits, Flathead Screwdriver, Dental Pick, Spectra Scissors, Needle Nose Pliers & Reel Lubricant. Extra space is allocated for a Click® Reels and a lighter, which can be purchased separately.	

RevoFit®

SKU	Product	Description	Picture
RF-200-07-01	RevoFit® Lamination Kit with Click® Reel	This single RevoFit Lamination Kit with the Click Reel allows you to fabricate user-controlled adjustability anywhere on a device. The Click Reel innovation creates adjustable devices with a lower streamlined profile, micro/macro adjustment that incrementally tightens and loosens with increased durability. Kit includes: 1x Click Reel, 1x Lamination Collar, 2m HD Lace, 1x Metal Lace Feeder, 3x Plastic Lace Feeders, 1x Reel Tool, 1.8m Lamination Tubing User Weight Limit: 300lb (135kg)	

RF-200-07-10	RevoFit® Lamination Fabricators 10-Pack with Click® Reel	The RevoFit Lamination Fabricators 10-Pack with the Click Reel is recommended for central fabricators and practices that create large volumes of adjustable devices. Save 10% on purchase price of each single kit while assuring you will always have a kit on hand.	
RP-002	HD Lace Kit	2m of HD (High Durability) Lace and Click Reel Components for RevoFit Lamination Kits with the Click Reel.	
RP-002-07	Click® Reel and HD Lace Kit	Replacement Click Reel and 2m of HD (High Durability) Lace for RevoFit Kits only. Note: This is not for replacing a BOA® dial. Please see RP0100-320-05 for replacing a BOA® dial.	
RP-300	Slider Part Kit	Replacement Slider and Lock Wire for dynamic AFOs, large gap prostheses or RevoLock applications.	

RP-309	Slider Lamination Kit	Lamination Kit for fabricating new dynamic AFOs, large gap prostheses or replacement for RevoLock Kits.	
RP-200	Lamination Collar	Replacement Lamination Collar for RevoFit & RevoLock Lamination Kits with the Click Reel only.	
RP-201	Lamination Tube - 1.8m	Replacement Lamination Tube (1.8m) for all Revo Kits with Click Reel and BOA® dials.	
RP-300	Slider Part Kit	Replacement Slider and Lock Wire for dynamic AFOs, large gap prostheses or RevoLock applications.	
RP-309	Slider Lamination Kit	Lamination Kit for fabricating new dynamic AFOs, large gap prostheses or replacement for RevoLock Kits.	
BL-002	HD Lace - 25m Spool	25m spool of HD (High Durability) Lace assures you will always have replacement lace on hand. Cost effective and sustainable with less waste. Compatible with both the RevoFit Lamination Kit with Click Reel and RevoFit Lamination Kit with BOA® dial.	
AC-100-05	Comfort Grip Ring - Black 5-Pack	Provide a soft-touch feel and 10% increased turning power to the existing Click Reel. Great for patients with dexterity issues.	

AC-101-05	Power Grip Ring - Black 5-Pack	Provide a soft-touch feel and 20% increased turning power to the existing Click Reel. Great for patients with dexterity and hand strength issues.	
AC-000-00	Click Care Kit™	The Click Care Kit enables you to service your clients in the office or on the road with everything readily accessible for maintaining, repairing or replacing Click components. Contents: 1x Surface Collar Kit, 2x HD Lace Kits, 2x HC Lace Kits, Flathead Screwdriver, Dental Pick, Spectra Scissors, Needle Nose Pliers & Reel Lubricant. Extra space is allocated for a Click® Reels and a lighter, which can be purchased separately.	

RevoLock®			
SKU	Product	Description	Picture

RL-301-07-01	RevoLock® Align Kit with Click® Reel	This single RevoLock Align Kit with the Click Reel is for transferable alignment fabrication for both upper and lower extremity sockets. The Click Reel innovation offers increased durability and power, a lower streamlined profile, and micro/macro adjustment that incrementally tightens and loosens. The Slider enhances safety for limbs changing in position or fluctuating in volume. Kit includes: 1x Click Reel, 1x Lamination Collar, 1x Slider Lamination Kit, 2x RevoLock LE Quick Connect Inserts, 1x RevoLock LE Quick Connect Snap, 2m HC Lace, 1x Metal Lace Feeder, 3x Plastic Lace Feeders, 1x Reel Tool, 0.5m Lamination Tubing, 1x Align Housing + Fabrication Post User Weight Limit: 300lb (135kg)	 The image displays the components of the RevoLock Align Kit. It includes a blue and black box labeled 'RevoLock ALIGN KIT'. The components laid out are: a black lamination collar, a black click reel, a black slider lamination kit, two black RevoLock LE Quick Connect Inserts, one black RevoLock LE Quick Connect Snap, a coiled green HC lace, one metal lace feeder, three plastic lace feeders, one black reel tool, a coiled black lamination tubing, and a black align housing with a fabrication post.
RL-300-07-01	RevoLock® 4-Hole Kit with Click® Reel	This single RevoLock 4-Hole Kit with the Click Reel is for standard alignment fabrication for lower extremity sockets and comes with a standard 25mm x 65mm distal housing and integrated 4-Hole plate. The Click Reel innovation offers increased durability and power, a lower streamlined profile and micro/macro adjustment that incrementally tightens and loosens. The Slider enhances safety for limbs changing in position or fluctuating in volume. RevoLock® Tool Form required for fabrication of RevoLock 4-Hole kit. Kit includes: 1x Click Reel, 1x Lamination Collar, 1x Slider Lamination Kit, 2x RevoLock LE Quick Connect Inserts, 1x RevoLock LE Quick Connect Snap, 2m HC Lace, 1x Metal Lace Feeder, 3x Plastic Lace Feeders, 1x Reel Tool, 0.5m Lamination Tubing, 1x 4-Hole Puck, 1x Tool Form Plug User Weight Limit: 300lb (135kg)	 The image displays the components of the RevoLock 4-Hole Kit. It includes a blue and black box labeled 'RevoLock 4-HOLE KIT'. The components laid out are: a black lamination collar, a black click reel, a black slider lamination kit, two black RevoLock LE Quick Connect Inserts, one black RevoLock LE Quick Connect Snap, a coiled green HC lace, one metal lace feeder, three plastic lace feeders, one black reel tool, a coiled black lamination tubing, a black 4-hole puck, and a small black tool form plug.

RL-310-00-01	RevoLock® Tool Form	The 25mm x 65mm RevoLock Fabrication Tool Form is a required, reusable fabrication tool for RevoLock 4-Hole Kits.	
RL-302-22-01	RevoLock Lanyard V3 with BOA® Dial	This kit offers the push to engage and pull to release the interface of the BOA® Dial in a RevoLock Lanyard solution. It includes a shorter and smaller diameter distal housing which accepts the Quick Connect Snap for easy connection to the liner. Kit includes: 1x BOA® Lanyard Dial, 1x Lamination Collar, 1x 0.5m Lamination Tubing, 2x LE Quick Connect Inserts, 1x Quick Connect Snap, 1x Align Housing + Fabrication Post, 2x 1.22m HC Lace, 1x Metal Lace Feeder User Weight Limit: 220lb (100kg)	
RP-310	RevoLock V2 & V3 BOA® Dial and HC Lace Kit (Port C)	Replacement 1x BOA® Lanyard Dial, 1.22m of HC Lace (High Capacity), 1x Metal Lace Feeder. For use with Port C Lace exit. Compatible only with RevoLock Lanyard V2 and V3 Kits	
RP-310A	RevoLock V2 & V3 BOA® Dial and HC Lace Kit (Port A)	Replacement 1x BOA® Lanyard Dial, 1.22m of HC Lace (High Capacity), 1x Metal Lace Feeder. For use with Port A Lace exit. Compatible only with RevoLock Lanyard V2 and V3 Kits	

RP-312	RevoLock V2 & V3 HC Lace Kit	1.22m of HC Lace (High Capacity), 1x Metal Lace Feeder. Compatible only with RevoLock Lanyard V2 and V3 Kits	
RP-001	HC Lace Kit	2m of HC (High Capacity) Lace and Click Reel components for RevoLock Align, RevoLock 4-Hole, and RevoSurface Tool Kits with the Click Reel.	
RP-000-07	Click® Reel and HC Lace Kit	Replacement Click Reel and 2m of HC (High Capacity) Lace for RevoLock Align, RevoLock 4-Hole and RevoSurface Tool Kits only. This is not for replacing a BOA® dial. Please see RP0100-320-05 for replacing a BOA® dial.	
RP-200	Lamination Collar	Replacement Lamination Collar for RevoLock Align, RevoLock 4-Hole, RevoLock Lanyard V3 & RevoFit Lamination Kits.	

RP-201	Lamination Tube - 1.8m	Replacement Lamination Tube (1.8m) for all Revo Kits with Click Reel and BOA® dials.	
RP-302	RevoLock® LE Quick Connect Insert	Replacement RevoLock LE Quick Connect Insert for connection to lower extremity liners. Compatible with RevoLock LE Quick Connect Snap (RP-303) found in the RevoLock Align, RevoLock 4-Hole and RevoLock Lanyard V3 Kits. Often purchased for outfitting additional patient liners.	
RP-303	RevoLock® LE Quick Connect Snap	Replacement RevoLock LE Quick Connect Snap connects the HC (High Capacity) Lace to the liner via the RevoLock LE Quick Connect Insert (RP-302).	
RP-305	RevoLock® Tool Form Plug 5-Pack	The RevoLock Tool Form Plug 5-Pack is to ensure fabricators have extra plugs on hand.	
RP-306	RevoLock® UE Fixed Insert	RevoLock Fixed Threaded Insert (M6) for upper extremity and pediatric liners. Compatible with RevoLock 4-Hole, RevoLock Align and RevoLock Lanyard V3 Kits. Best for patients who cannot operate the RevoLock Quick Connect Snap or for suction systems needing an airtight seal. This part is NOT releasable and does not connect to the RevoLock LE Quick Connect Snap.	

RP-307	RevoLock® LE Fixed Insert	RevoLock Fixed Threaded Insert (M10) for lower extremity liners only. Compatible with RevoLock 4-Hole, RevoLock Align and RevoLock Lanyard V3 Kits. Best for patients who cannot operate the RevoLock Quick Connect Snap or suction systems needing an airtight seal. This part is NOT releasable and does not connect to the RevoLock LE Quick Connect Snap.	
RP-300	Slider Part Kit	Replacement Slider and Lock Wire for dynamic AFOs, large gap prostheses or RevoLock applications.	
RP-309	Slider Lamination Kit	Lamination Kit for fabricating new dynamic AFOs, large gap prostheses or replacement for RevoLock Kits.	
BL-001	HC Lace - 25m Spool	25m spool of HC (High Capacity) Lace assures you will always have replacement lace on hand. Cost effective, sustainable and less waste. Used in the RevoSurface, RevoLock Align, RevoLock 4-Hole, and RevoLock Lanyard V2 with BOA Dial.	
AC-100-05	Comfort Grip Ring - Black 5-Pack	Provide a soft-touch feel and 10% increased turning power to the existing Click Reel. Great for patients with dexterity issues.	
AC-101-05	Power Grip Ring - Black 5-Pack	Provide a soft-touch feel and 20% increased turning power to the existing Click Reel. Great for patients with dexterity and hand strength issues.	


AC-000-00

Click Care Kit™

The Click Care Kit enables you to service your clients in the office or on the road with everything readily accessible for maintaining, repairing or replacing Click components.

Contents: 1x Surface Collar Kit, 2x HD Lace Kits, 2x HC Lace Kits, Flathead Screwdriver, Dental Pick, Spectra Scissors, Needle Nose Pliers & Reel Lubricant. Extra space is allocated for a Click® Reels and a lighter, which can be purchased separately.





QuickFit® (Made with the BOA® dial)

SKU	Product	Description	Picture
OK0180-130-05	QuickFit® Lacer, XS	The extra small size QuickFit Lacer is designed to close a small region of an orthosis, a pediatric sized orthosis, or pediatric footwear up to size 35 EU. Comes pre assembled with a low-profile dial and 4 guides. QuickFit Rivets sold separately, but are required to install the QuickFit Lacer kits. Lacers are sold individually. For a pair of shoes, two kits are required.	
OK0185-130-05	QuickFit® Lacer, XS 10-Pack	The extra small size QuickFit Lacer is designed to close a small region of an orthosis, a pediatric sized orthosis, or pediatric footwear up to size 35 EU. Comes pre assembled with 10 extra small QuickFit Lacer kits and 100 QuickFit Rivets. For a pair of shoes, two kits are required.	N/A
OK0200-130-05	QuickFit® Lacer, SM	The small size QuickFit Lacer is designed to close a small region of an orthosis, a small sized orthosis, or footwear size 36 - 41 EU. Comes pre assembled with a low profile dial and 6 guides. QuickFit Rivets sold separately, but are required to install the QuickFit Lacer kits. Lacers are sold individually. For a pair of shoes, two kits are required.	
OK0205-130-05	QuickFit® Lacer, SM 10-Pack	The small size QuickFit Lacer is designed to close a small region of an orthosis, a small sized orthosis, or footwear size 36 - 41 EU. It comes with 10 Small QuickFit Lacer kits and 100 QuickFit Rivets. For a pair of shoes, two kits are required.	N/A

OK0220-130-05	QuickFit® Lacer, MD	The medium size QuickFit Lacer is designed to close a medium sized orthosis or footwear size 42 EU and larger. Comes pre assembled with a low profile dial and 8 guides. QuickFit Rivets sold separately, but are required to install the QuickFit Lacer kits. Lacers are sold individually. For a pair of shoes, two kits are required.	
OK0225-130-05	QuickFit® Lacer, MD 10-Pack	The medium size QuickFit Lacer is designed to close a medium sized orthosis or footwear size 42 EU and larger. It comes with 10 Medium QuickFit Lacer kits and 100 QuickFit Rivets. For a pair of shoes, two kits are required.	N/A
OK0240-230-05	QuickFit® Lacer, LG	Designed to close or secure above the ankle shoes or orthoses, such as a Leather AFO's, plastic rigid AFO's, or Crow boots. Can also be used on footwear such as hiking boots. Comes pre assembled with a more powerful dial and 10 guides. Rivets sold separately.	
OK0165-230	QuickFit® Releasing Lacer, LG	Adds opening functionality to our QuickFit Lacer by adding releasing guides to one side so that the system can fully open like a door for easy donning. Comes with a more powerful reel to close larger devices. This kit is unassembled with 4 fixed guides and 3 releasing guides. Rivets sold separately.	

TL-001	QuickFit® Tool	The QuickFit Tool is customized for use with QuickFit Rivets and is the recommended installation tool. The QuickFit Tool properly secures rivets and provides either an 1/8" and 3/16" rivet hole.	
SP1040	QuickFit® Rivets 100-Pack	QuickFit Rivets come in packs of 100. You can supply your own rivets, or purchase our easy to use rivet pack to install QuickFit Lacers to any device or shoe.	
RP0200-130-05	QuickFit® XS/SM/MD Dial and Lace Replacement Kit	Replacement BOA® dial and 1.5m of Low-Mid Power Steel Lace for XS/SM/MD size QuickFit Lacers using BOA® dials.	
RP0200-230-05	QuickFit® LG Dial and Lace Replacement Kit	Replacement BOA® dial and 2m of compatible Low-Mid Power Steel Lace for LG size QuickFit Lacers using BOA® dials.	
LOK019	Low-Mid Power Steel Lace - 10m Spool	Bulk replacement Low-Mid Power Steel Lace (10m spool) compatible only with BOA® dials and QuickFit Lacer kits. NOT COMPATIBLE with Click® Reel.	

OK1410-130-05	QuickFit® Strap 1.5"	QuickFit Straps are designed to quickly install on a wide variety of cuffs, orthoses, and prosthetic devices to create an improved, easy to use, adjustable strap. • Prefabricated, durable, and installs in minutes • Don/doff via simple Velcro® closure • 76mm (3") of adjustable travel • Usable strap length: Max: 0.75m (30") / Min: 0.2m (9") • Cut to size • Width (38mm)	
OK1415-000-05	QuickFit® Buckle 1.5"	The QuickFit Buckle can be added onto the QuickFit Strap in a few seconds. The magnetic buckle allows the user to easily connect the strap - even with limited dexterity or in hard to reach areas. Its plastic on plastic interface eliminates wear issues of Velcro® closure. • Easily adds onto QuickFit Strap in seconds • Adds 76mm (3") to strap usable length • 5mm (.2") profile • Width (38mm)	

Legacy Replacement Parts using the BOA® dial

SKU	Product	Description	Picture
PK0060-000-10	RevoFit® Re-Lacing Kit for BOA® dial	2m of HD (High Durability) Lace for RevoFit Lamination Kits Compatible only with RevoFit Kits with BOA® dial.	
RP0100-320-05	RevoFit® Dial and Lace Replacement Kit with BOA® dial	Replacement BOA® dial and 2m of HD (High Durability) Lace Compatible only with RevoFit Kits with BOA® dial.	
RP3000-220-05	RevoLock® Lanyard Dial and Lace Replacement Kit with BOA® dial	Replacement BOA® dial and 1.2m of lace Compatible only with RevoLock Kits with BOA® dial. Not compatible with RevoLock V2 and V3 Kits.	

RP3001-000-05	RevoLock® Lanyard Re-Lacing Kit for BOA® dial	Contains 1.2m of lace Compatible only with RevoLock Kits with BOA® dial. Not compatible with RevoLock V2 and V3 Kits.	
RP3002-000-05	RevoLock® Lanyard - Lower-Ex Threaded Insert (10mm)	Replacement RevoLock Lower-Extremity Threaded Insert which attaches the lace to the liner via a threaded connection. Compatible only with RevoLock Lanyard & RevoLock Lanyard V2 Kits with BOA® dial.	
RP3003-000-05	RevoLock® Lanyard - Upper-Ex Threaded Insert (6mm)	Replacement RevoLock Upper-Extremity Threaded Insert which attaches the lace to the liner via a threaded connection. Compatible only with RevoLock Lanyard & RevoLock Lanyard V2 Kits with BOA® dial.	