

New QuickFit® Lacer

Available July 1, 2026



The Most Versatile Dial-Based Closure System

Add One-Handed Adjustability to Dial in the Perfect Fit

WHAT'S NEW FOR 2026

- ✓ **New kit design** offers greater customization and flexibility
- ✓ **NEW dual-function lace guides** expand routing possibilities
- ✓ QuickFit Lacer now available in **White**



Partially Pre-Assembled

Lace pre-installed in one port of the dial makes custom lace configuration easier.



Dual-Function QuickLace Guides

Operate in fixed or releasable modes for flexible configurations and easier donning and doffing.



Steel Lace

Nylon-coated steel lace is low friction and high strength.



Easily Replaceable Components

Dial and lace replacement is easy, with no need to remove rivets or re-fabricate device.

QuickFit® Lacer: A low-profile dial system designed for performance and comfort with smooth, precise adjustments

The QuickFit Lacer is a dial-based lacing system that tightens a high-strength lace to provide adjustable closure for upper extremity prostheses, orthoses, and footwear. The low-profile dial enables smooth, one-handed micro-adjustments while distributing tension evenly across the closure area.

Redesigned for easy integration, the system comes partially pre-assembled and allows practitioners to customize lace routing by placing guides where needed—no system disassembly required. QuickFit Lacers are sold in pairs, making them ideal for footwear and bilateral applications.

Available in black or white. Additional QuickLace Guides are available separately.

QuickFit® Lacer Black 2-Pack

QFL-130-B-2



QuickFit® Lacer White 2-Pack

QFL-130-W-2



QuickFit® Lacer Large: Designed for above-the-ankle footwear and orthoses

The QuickFit Lacer Large is a dial-based lacing system that tightens a high-strength lace to deliver powerful, adjustable closure for above-the-ankle footwear and orthoses. Designed specifically for applications such as leather AFOs, rigid plastic AFOs, and CROW boots, it provides dependable closure with a simple turn of a dial.

Redesigned for easy integration, the system comes partially pre-assembled and allows practitioners to customize lace routing by placing guides where needed. QuickFit Lacers are sold in pairs, making them ideal for bilateral applications.

Available in Black. Additional QuickLace Guides are available separately.

QuickFit® Lacer Large Black 2-Pack

QFL-230-B-2



QuickFit® Lacer Product Parts



QuickFit® Tool

TL-001



QuickFit® Rivets Silver, 100-Pack

SP1040-S



QuickFit® Rivets Black, 100-Pack

SP1040



QuickFit® Dial and Lace Replacement Kit - White

RP0200-130-W



QuickFit® Dial and Lace Replacement Kit - Black

RP0200-130-05



QuickFit® LG Dial and Lace Replacement Kit - Black

RP0200-230-05



Steel Lace - White 10m Spool

LOK019-W



Steel Lace - Black 10m Spool

LOK019



QuickLace Guide Black 10-Pack

RP-103-B



QuickLace Guide White 10-Pack

RP-103-W

SCAN FOR INSTRUCTIONS

www.clickmedical.co/instructions



For additional tech specs, please go to:
www.clickmedical.co/quickfit-lacer